Japet.^{w+}

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Table of contents

About this guide	30
Properties and copyright	31
01. Introduction and warnings	32
Presentation of the device	32
Clinical performance of the device	33
Indications	33
Contraindications	33
Vigilance	33
Warning	33
Contents of the case	34
	34
02. Using Japet.W⁺	
Measurements	36
How the system works	36
Battery and battery installation	37
Installation: second skin and active system	37
Tightening	38
Starting up	40
Removal and folding	41
Cleaning	42
Battery charging	43
Disposal of the device and components	44
	44
03. Technical characteristics	
Device features	45
Battery electrical features	45
Charger electrical features	45
Mechanical features	45
Performance claims	46
Connectivity features	46
	47

04. Conditions of use, storage and transport	46
Operating conditions	46
Storage and transport conditions	46
Electromagnetic compatibility	47
05. Troubleshooting	48
06. Regulatory information	48
Liability	48
Guarantees	48
Trademarks	49
Patents	49
07. Other information	50
Symbols on the label and packaging	50
Version history of the user guide	51
Contact information	52
Contacts	52

About this guide

This user's guide has no contractual value and JAPET MEDICAL DEVICES® cannot be held responsible in any way for the information contained in this guide.

This user's guide details all the knowledge required to implement, use and maintain the JAPET.W device+ and the collection of information displayed.

Thus, after a careful and complete reading before using the device, the operator will be able to:

- Connect the peripheral elements (power cord) and turn on the JAPET.W device+.
- Configure the device, perform basic maintenance.

JAPET MEDICAL DEVICES® provides this guide "as it is", without warranty of any kind, either expressed or implied, including, but not limited to, the implied warranties of merchantability and/or fitness for a particular purpose in order to provide simple and accurate information. JAPET MEDICAL DEVICES® therefore assumes no responsibility for any misinterpretation. Although every effort has been made to ensure that the manual is as accurate as possible, it may contain technical inaccuracies and/or typographical errors. In no event shall JAPET MEDICAL DEVICES® be liable for any loss of profit, loss of business, loss of data, business interruption, or for indirect, special, incidental, or consequential damages of any kind.

In the event of damage arising from a defect (imperfection) or error contained in this user guide, JAPET MEDICAL DEVICES® undertakes to send the customer a hard copy or electronic document containing the corrections made to this guide as soon as possible.

This guide is regularly updated. The most recent version of this guide is available on request from JAPET MEDICAL DEVICES®. However, in the event that major changes are made to the guide, JAPET MEDICAL DEVICES® undertakes to send the customer the new guide in paper or electronic form as soon as possible. It is specified that this does not imply the updating of the hardware and/or software in your possession.

It is the responsibility of the owner of the product to keep this manual for as long as the product is in use. This guide includes a chapter on a quick solution to the most frequently encountered problems. Any request for information or modification of this guide should be addressed to

**JAPET MEDICAL DEVICES ®,
F3A
12, avenue Pierre Mauroy
59120 LOOS
France**

Properties and copyright

All manuals and documentation of any kind are the property of JAPET MEDICAL DEVICES® and are protected by copyright, all rights reserved. Your right to copy this documentation is limited to the rights of legal copying. These manuals may not be distributed, translated or reproduced, in whole or in part, in any manner or form, without the prior written permission of JAPET MEDICAL DEVICES®. Thus, reproduction, adaptation or translation of this guide without prior written permission is prohibited to the extent permitted by copyright laws.

Copyright© - 06/2018 - JAPET MEDICAL DEVICES® - All rights reserved

01. Presentation and warnings

Presentation of the device

Low back pain is a major public health problem that is constantly increasing in Western countries where it affects or will affect nearly 4 out of 5 people.


Among its main causes, the crushing of intervertebral discs can cause anything from a simple persistent pain to paralysis. This pain, when it persists and intensifies, creates disabling pain that restricts physical activities and gradually isolates the patient from his social and professional environment. The complexity of the origin of the pathology and the multidisciplinary management amplifies the risks of a transition to chronicity. JAPET.W⁺, the exoskeleton to relieve the pain of low back pain, while allowing full mobility.

This is based on devices commonly used in rheumatology such as elongation tables to relieve the patient by decompressing the spine.


It is produced by the thrust of micromotors supported on two belts, one on the pelvis, the other on the chest.

The device is also inspired by robotic exoskeletons in order to adjust the position of the micromotors according to the user's movements to free up trunk mobility. This allows the user to resume the physical activity recommended for this type of pathology.

Active system

Second skin

Clinical performance of the device

Decompression of the intervertebral discs.

Indications

The JAPET.W⁺ device is indicated for people with chronic and acute low back pain. The device is intended for ambulatory distraction of the lumbar spine.

It is recommended to use the device over a T-shirt.

The use of the device is intended for healthcare professionals and/or referents who have been trained by JAPET MEDICAL DEVICES® or by a third party authorised to carry out this training.

It is recommended to use the device in 30-60 minute sessions (2-3 times a day).

For use in the workplace, it is recommended that the JAPET.W⁺ device be worn before a painful episode.

Contraindications

This device is contraindicated in the following cases:

- Neurological root and spinal cord deficits
- Treatment with implanted neurostimulation
- Heart or circulatory disease and severe respiratory problems
- Tumour or infectious process of the vertebrae
- Osteoporosis (< -2.5 ds)
- History of arthrodesis surgery
- Herniated disc surgery (< 3 months)
- Fracture of the dorsolumbar spine (< 3 months)
- Rib fracture (< 3 months)
- Skin injury, bruising and stretching injuries of the trunk
- Pregnancy (> 4 months)
- Previous lumbar prosthesis surgery

If any of the above conditions apply to you, consult your doctor before using the device.

Possible adverse effects of the device

Incorrect positioning of the device or incorrect sizing can cause :

- Haematoma,
- Irritation,
- Pain.

In case of an adverse reaction not listed and related to the device, please contact JAPET MEDICAL DEVICES®.

Contact information is available at the end of the user guide.

Vigilance

In the event of a serious incident, i.e. any incident that directly or indirectly resulted in, may have resulted in or may result in :

- The death of a patient, user or any other person;
- A serious temporary or permanent deterioration in the state of health of a patient, user or any other person;
- A serious threat to public health;

Please contact JAPET MEDICAL DEVICES® and the relevant authorities.

Contact information is available at the end of the user guide.

Warnings


-  Do not use other batteries than those distributed by JAPET MEDICAL DEVICES®.
-  Do not use other battery chargers than the one distributed by JAPET MEDICAL DEVICES®.
-  To avoid the risk of electric shock, the battery charger should only be connected to a mains supply with a protective earth.
-  If the protective case is damaged, please contact JAPET MEDICAL DEVICES®. Do not use the JAPET.W⁺ device if it is damaged in any way.
-  Take routine precautions to avoid contact with blood or other bio-contaminating materials. Contaminated material should be treated in accordance with the hygiene requirements of the establishment.
-  Do not expose the device to rain, dripping or splashing liquids. Do not immerse it in water.
-  The device is internally powered by a battery which is recharged externally with the supplied charger.
-  It is strictly forbidden for a third party, other than a technician authorised by JAPET MEDICAL DEVICES®, to open or modify the device.
-  It is strictly forbidden to dismantle or modify the device (this act interrupts all guarantees).
-  The motors may become hot during normal use of the device, the maximum temperature that can be in contact with the skin is 41°C (105.8°F). In case of abnormal or excessive heating of the motors, stop wearing the device and contact JAPET MEDICAL DEVICES®.

Contents of the case

- 1 Active system
- 2 Second Skins
- 1 User Guide
- 2 RRC 2040 batteries
- 1 Suitable charger and its power cable
- 2 Wash bags

Active system


Front view

Right side view


Back view


Second skin


External part


Inner part


Accessories


CHARGER

WASHING BAG (x2)

BATTERY (x2)

USER GUIDE

02. Using Japet.W+

Measurements

There are two different sizes of active system (T1 and T2), as well as four sizes of second skin (S, M, L, XL).


In order to select the correct sizes, measure the underbust and hip circumferences as shown below.

Now plot the values in the following graph


to find the size of the device that fits the user's body type.

The device must be used with a second skin adapted to the active system and to the user's morphology.

Schemas

Graph

Second skin

S	UCC	80 - 89 max
	HC	86 - 96 max
M	UCC	89 - 98 max
	HC	96 - 106 max
L	UCC	98 - 110 max
	HC	106 - 120 max
XL	UCC	110 - 122 max
	HC	120 - 134 max

Active system

	Size 1
	Size 2



How the system works

Follow these steps consecutively to use the JAPET.W⁺ device.

Battery and battery installation

a. Check the state of charge of the battery with the indicator light.

b. Install the battery in the case, with the flap and indicators facing outwards, until you hear a beep.

Battery level indicators

	Battery full	Four LEDs light up gradually, the battery level is 100%.
	Battery almost full	Three LEDs light up gradually, battery level is 75%.
	Battery half exhausted	Two LEDs light up gradually, battery level is 50%.
	Battery almost exhausted	Only one LED light up, battery level is 25%.
	Battery exhausted	One LED flashes, the battery level is below 5%. Recharge the battery according to the instruction.

Installation: second skin and active system


There are two different ways to connect the second skin to the active system:

- option 1 : on a table
- option 2 : on the body


Preferably option 1.

Option 1: on a table.

a. Place the active system on a flat surface with the label visible. Then attach the second skin to the active system by overlapping the orange velcro. Ensure that all parts of the second skin and the active system are overlapping.

b. Position the device around your waist and close the fasteners. Check that the device is placed symmetrically on the hips. To do this, check the alignment of the lines on the upper and lower belts. If necessary, reposition the device.

Option 2: on the body

a. Put the second skin around the waist so that the orange velcro is on the outside, the label on the back (inside, top).


b. First fasten the top bindings together and then the bottom ones. Do not overtighten the second skin, the active system rollers will allow it to be adjusted to the body shape.

c. Check that the seam line of the trousers and the marker line of the active system are aligned so that the assembly is properly positioned on the hips. Readjust if necessary.

d. Take out the active system and then hold the end of the two fasteners where the clips are located in your hand. Slip it over the second skin so that the side with the orange velcro fasteners is facing inwards.


To properly attach the second skin part to the active system part, make sure that the orange velcro overlaps.

e. Check that the device is placed symmetrically on the hips. To do this, check the alignment of the lines on the upper and lower belts. If necessary, reposition the device.


Tightening

a. Press in the reel (until you hear a click).
Once pressed, turn the reel clockwise.


b. Repeat this process for the top and bottom until a comfortable fit is achieved.

Switching on

1. Press the ON/OFF button for 1 second.
2. To increase the traction level press the + button.

3. To decrease the traction level, press the - button. An audible signal confirms the pressure. Check that the device does not slip when deploying the actuators.




Removal and folding

a. Once the traction level is at 0, pull on the reels to release the device and press the ON/OFF button.


b. Remove the active system from the second skin and place it in its storage box, holding the micromotors towards you.

c. Then remove the second skin. Fold it in 4 from left to right and then from top to bottom and store it in its wash bag and then its storage case.

Cleaning



We recommend to clean the second skin at 40° maximum in the washing bag present in the suitcase, after a maximum use of 2 weeks. Drying can be done in the open air.

Battery charging

When the battery level is insufficient (about 5% left), two indicators are activated :

- A repetitive beep
- An indicator light (flashing) at the first white LED (fig.a)
- Then remove the battery from the device and place it in the charger as shown. (fig.b)

Ensure that the charger is connected to the mains using the cable provided.

- A red light indicates that the battery is charging.
- A green light comes on when charging is complete.

 **WARNING:** To avoid the risk of electric shock, the battery charger should only be connected to a mains supply with a protective earth.

Disposal of the device and components

The Japet.W⁺ device is subject to the following directives:

- Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment.
- Directive 2006/66/EC of the European Parliament and of the Council of 6 September 2006 on batteries and accumulators and waste batteries and accumulators

The Japet.W⁺ device contains (lithium) batteries, electrical wires and electronic components that cannot be disposed of with household waste.

Do not dispose of the batteries supplied by JAPET MEDICAL DEVICES® or the device with normal household waste at the end of the product's life, but take them to an appropriate recycling centre, thus helping to preserve the environment.

03. Technical specifications

Device features

Manufacturer	JAPET MEDICAL DEVICES® F3A – 12, avenue Pierre Mauroy 59120 LOOS FRANCE
Model	JAPET.W ⁺
Medical device classification	Class IIa according to regulation 9, medical device regulation
Type of party applied	BF (the whole JAPET.W ⁺ device)
IP code	IP21
Embedded software class	Class A
Lifetime	5 years

Battery electrical features

Model	RRC2040
Rated voltage	10,8V
Capacity	3350mAh
Dimension	84,9mm x 58,8mm x 21,9mm
Mass	170 g
Lifetime at 25°C	> 300 cycles with a minimum of 75% of its initial capacity

Charger electrical features

Model	RRC-SMB-MBC
Input voltage	100 - 240VAC
Input current	2,80A max
Consumption at rest	No load <0.5 W @ 230 VAC

Mechanical properties

Dimensions	104cm x 30cm x 10cm
Mass of the device	<2,4 Kg

Performance claims

Actuator thrust	4Kg
Overall thrust	16Kg

Connectivity features

Connectivity	Bluetooth BLE
Receiving frequency band	2400 MHz to 2483.5 MHz
Transmission frequency band	2400 MHz to 2483,5 MHz
Type	BLE modulation
Equivalent radiated power	Value below 0 dBm

04. Conditions of use, storage and transport

Operating conditions

Temperature	+0°C to +40°C (+32°F à +104 °F)
Humidity in the air	5% to 95%

Storage and transport conditions

Temperature	-20°C to +40°C (-4°F to +104°F)
Humidity in the air	5% to 95%
Atmospheric pressure	57 kPa to 106 kPa

Electromagnetic compatibility

Medical electrical equipment is subject to special precautions regarding electromagnetic compatibility (EMC) and should only be installed and used in accordance with the information specified in this document.

Wireless communication devices such as wireless home network routers, mobile phones, cordless phones and their cradles, and walkie-talkies may interfere with JAPET.W⁺. Other cables and accessories not supplied by the manufacturer may impair the EMC performance and cause increased emissions or decreased electromagnetic immunity of JAPET.W⁺.

The JAPET.W⁺ device should be installed

 **WARNING:** Use of this unit next to or stacked with other units should be avoided because it may cause malfunction. If such use is necessary, observe this unit and other units for normal operation.

 **WARNING:** The use of accessories, transducers, and cables other than those specified or provided with this device may result in increased electromagnetic emissions or decreased immunity of this device and may result in improper operation.

and operated in hospitals and away from HF sources of magnetic resonance imaging. The performance of the JAPET.W⁺ device corresponds to an overall thrust of 16 Kg. In case of electromagnetic disturbances this performance could be lost and the therapeutic effectiveness of the device could be limited.

The JAPET.W⁺ device complies with the IEC 60601-1-2 : 2014 (Ed.4) standard.

 **WARNING:** Portable RF communications devices (including peripherals such as antenna cables and external antennas) should not be operated closer than 30 cm (12 inches) to any part of the device. Otherwise, the performance of this device may be impaired.

05. Troubleshooting

In case of malfunction, remove the device and remove the battery. Make sure the battery is charged and then put it back into the device. Press the "⏻" button to restart the device. If the malfunction persists,

contact JAPET MEDICAL DEVICES®.

Do not disassemble the device.

If you need new batteries, please contact JAPET MEDICAL DEVICES®.

06. Regulatory information

Responsibility

JAPET MEDICAL DEVICES® cannot be held responsible for any damage that may result from the use of this product under conditions that do not conform to the instructions given in this guide.

MEDICAL DEVICES® shall authorise the purchaser to return all or part of the guaranteed product.

Guarantees

JAPET MEDICAL DEVICES® warrants to the purchaser that the device ("the warranted product") is free from defects in workmanship and materials under normal, proper and intended use.

JAPET MEDICAL DEVICES® will decide whether to repair or replace the product and parts covered by this warranty, and all replaced products or parts will become the property of JAPET MEDICAL DEVICES®.

The obligations of JAPET MEDICAL DEVICES® under this warranty are to repair or replace all or part of the warranted product which JAPET MEDICAL DEVICES® deems reasonable to cover as defective, provided that the purchaser asserts the warranty during the warranty period and the product is returned to JAPET MEDICAL DEVICES® in accordance with the provisions of the general conditions of sale or distribution agreements.

During the warranty period, JAPET MEDICAL DEVICES® may, but shall not be obliged to, make engineering improvements to all or part of the warranted product. If JAPET MEDICAL DEVICES® deems reasonable that a repair or replacement is covered by the warranty, JAPET MEDICAL DEVICES® will pay the cost of shipping the repaired or replacement product to the purchaser. Under this warranty, the risk of loss or damage during shipment shall be borne by the party shipping the product. Products returned by the purchaser under this warranty will be packaged in the original packaging or equivalent packaging to protect the product. If the buyer returns a product to JAPET MEDICAL DEVICES® in inadequate packaging, any damage to the material found on receipt by JAPET MEDICAL DEVICES® (and not already reported) shall be presumed to have occurred in transit and shall be the responsibility of the buyer.

Repair or replacement of products under this warranty does not extend the warranty period.

To request repair or replacement under this warranty, the purchaser should contact JAPET MEDICAL DEVICES® directly. JAPET

This warranty shall not extend to all or part of the guaranteed products which have been subject to misuse, negligence or accident; which have been damaged by causes external to the products, which have been used against the instructions of JAPET MEDICAL DEVICES®, which have been modified, dismantled, repaired or reassembled by a technician who is not a member of the JAPET MEDICAL DEVICES® staff. JAPET MEDICAL DEVICES® shall not be obliged to make any repairs, replacements or modifications resulting in whole or in part from normal wear and tear. JAPET MEDICAL DEVICES® does not grant its guarantee

- (a) any product that is not guaranteed
- (b) any product purchased from a supplier other than JAPET MEDICAL DEVICES® or its authorised distributor or;
- (c) any product sold under a brand name other than JAPET MEDICAL DEVICES®.

THIS WARRANTY IS THE SOLE AND EXCLUSIVE WARRANTY OF JAPET MEDICAL DEVICES® PRODUCTS MADE ONLY TO THE PURCHASER AND IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. JAPET MEDICAL DEVICES® LIABILITY ARISING OUT OF THE SALE OF THE PRODUCTS OR THEIR USE UNDER WARRANTY, CONTRACT, TORT OR OTHERWISE SHALL NOT EXCEED THE AMOUNTS ACTUALLY RECEIVED BY JAPET MEDICAL DEVICES®. JAPET MEDICAL DEVICES® SHALL NOT BE LIABLE FOR ANY DAMAGES, LOSSES OR EXPENSES OF A CASUAL, SPECIAL OR CONSEQUENTIAL NATURE (INCLUDING BUT NOT LIMITED TO LOST PROFITS) RESULTING DIRECTLY OR INDIRECTLY FROM THE SALE, FAILURE TO SELL, USE OR INABILITY TO USE ANY PRODUCT. EXCEPT AS OTHERWISE PROVIDED IN THIS WARRANTY, ALL PRODUCTS ARE PROVIDED ON AN "AS IS"

BASIS WITHOUT WARRANTY OF ANY KIND, EXPRESSED OR IMPLIED.

Trademarks

JAPET MEDICAL DEVICES® is a registered trademark of JAPET MEDICAL DEVICES®.

Patents

The JAPET.W+ device is covered by several patents in France and internationally.

07. Miscellaneous information

Symbols on the label and packaging :

	Indicates that the article is a medical device
	Symbol indicating potentially dangerous situations and actions
	The name and address of the manufacturer
	The date of manufacture of the device
	The serial number of the device
	The manufacturer's reference for the device
	Unique device identifier
	Consult the user guide carefully
	CE marking of the device
	DC power supply
10,8 V	Voltage of the power supply delivered by the battery
3A	Intensity of the power supply delivered by the battery
	Applied part type BF
	Symbol indicating the low and high temperature limit for storage and transport

	Fragile, handle with care
	Fears moisture
	Humidity limitation
	Shielded from the light
IP21	Protection class
	Separate waste collection system
	Very delicate wash at maximum temperature 40°C
	Money laundering prohibited
	Tumble drying not allowed
	No ironing

Version history of the user guide

Version	Date	Changes
1.0	10/2020	Creation based on the Atlas user guide.
2.0	11/2020	Correction of the history and a contraindication.
3.0	02/2021	Translation into German and addition of a contraindication. Updated regulatory references.
4.0	03/2022	Change of address, regulatory references, belt tightening device, technical specifications, after-sales contact.
5.0	07/2022	Clarification of the indication, addition of adverse effects, vigilance, clinical performance, modification/addition of symbols
6.0	08 /2022	Addition of an adverse reaction: Pain
7.0	04/2023	Update for Japet.W ⁺

07. Other informations

Contact information

	JAPET MEDICAL DEVICES® F3A 12, avenue Pierre Mauroy 59120 LOOS FRANCE	
-----------------------------------------------------------------------------------	-----------------------------------------------------------------------------------	-----------------------------------------------------------------------------------

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Copyright© – 06/2018 – JAPET MEDICAL DEVICES® – Tous droits réservés
Copyright© - 06/2018 - JAPET MEDICAL DEVICES® - Alle Rechte vorbehalten